

## Bind-It™

### Designed for Nuclear Medicine

#### The Next Generation of Radioactive Decontaminants

Reasons to use Bind-It Decontaminant Hand Soap every time you wash your hands:

- Will not dry your hands like ordinary lab soap
- Fresh, clean scent
- Economical & convenient
- Touch-free wall mount dispenser
- Meets the intent of ALARA
- Removes contamination every time you wash

When you use Bind-It Decontaminant Spray the difference is clear:

- **Will NOT cloud lead glass**
- Safely decontaminate Syringe Shields & L-Blocks
- Removes dirt, dust, oils and grease as well as Radioactive Contamination
- Used regularly it provides a protective barrier that aids in subsequent cleanup
- Very mild, pleasant scent



LTI047



LTI046



LTI030



LTI051

LTI051

Bind-It Decontamination Introductory Kit

- (1) Touch-free dispenser
- (1) Bind-It Hand Soap- 32 oz bottle
- (1) Bind-It Decontaminant Spray - 32 oz bottle

LTI046

Bind-It Decontaminant Hand Soap - **8 oz bottles**, Case of 12

LTI047

Bind-It Decontaminant Hand Soap - **32 oz bottles**, Case of 12

LTI030

Bind-It Decontaminant Spray Cleaner - **32 oz bottles**, Case of 12



LTI031

Decontaminant Wipes, 50/Package



LTIHCP

Bind-It Patient Care Packs allow patients to safely and effectively minimize and remove contamination from their isolation areas.

- The 8oz Hand Soap will accommodate many hand washings over the course of the patient's post treatment period and can be used in the shower
- The Ready-To-Use Spray is ready to use on the bathroom sink, bathtub/shower and floor
- The Concentrate is ideal for cleaning the sink, bathtub/shower and toilet bowl (the largest typical point of contamination)

Bind-It Patient Care Pack

- Includes:
- (2) 8oz bottle of Hand Soap
  - (1) 32oz bottle of Ready-To-Use Spray
  - (1) 8 oz bottle of Concentrate