


# Tungsten Shield for Medrad® Automated Solaris Injector



3mm Tungsten  
Height: 1 7/8"  
Width: 1 11/16"  
Inside Diameter: 1 1/4"  
Bottom Hole Diameter: 1/2"  
Weight: 0.7 Lbs

TS-MED Tungsten Shield for Medrad® Automated Solaris Injector with Screw



## Lead-Lined Storage Container

Useful for storing used syringes, alcohol wipes, etc., that may be contaminated with low-energy gamma or beta radiation residue prior to disposal

DJ-202 Lead-Lined Storage Container, .125" Lead

Exterior Dimensions: 7" H x 6" Dia

Interior Dimensions: 6.25" H x 5.5" Dia

DJ-202-25 Lead-Lined Storage Container, .25" Lead

Exterior Dimensions: 7" H x 6" Dia

Interior Dimensions: 5 7/8" H x 5 1/8" Dia

PC-201 Plastic Insert for DJ-202, 12/bx

