



Pinestar Technology, Inc.

Bind-It™

Designed for Nuclear Medicine

The Next Generation of Radioactive Decontaminants

- Will not dry your hands like ordinary lab soap
- Fresh, clean scent
- Touch-free wall mount dispenser
- Meets the intent of ALARA
- Removes contamination every time you wash



When you use Bind-It Decontaminant Spray the difference is clear:

- **Will NOT cloud lead glass**
- Safely decontaminate Syringe Shields & L-Blocks
- Removes dirt, dust, oils and grease as well as Radioactive Contamination
- Used regularly it provides a protective barrier that aids in subsequent cleanup
- Very mild, pleasant scent

- LTI051** Bind-It Decontamination Introductory Kit
(1) Touch-free dispenser
(1) Bind-It Hand Soap- 32 oz bottle
(1) Bind-It Decontaminant Spray - 32 oz bottle



- LTI046** Bind-It Decontaminant Hand Soap - **8 oz bottles**, Case of 12



- LTI047** Bind-It Decontaminant Hand Soap - **32 oz bottles**
Available in: Case of 12 Qty 6-11 Qty 1-5

- LTI030** Bind-It Decontaminant Spray Cleaner - **32 oz bottles**
Available in: Case of 12 Qty 6-11 Qty 1-5

Email: ptipa@pinestar.com • www.pinstar.com



Toll Free 1-800-682-2226 Phone 724-932-2121 Fax 724-932-3176